# KEMOT®

## Breathalyser MIE0214

## User Manual

English 1-10

Polski 11-22

Română 23-34

Deutsch 35-46

## Table of Contents

Introduction	2
Preparation for correct measurement	2
Description of parts	2
Power supply	3
Operation	3
Accessories	7
Error and troubles	8
Important notes	9
Specifications	9
Calibration and record	10
Precautions	10


## Introduction

The device is intended to measure alcohol in the human breath. Measurements obtained by this device are used in the diagnosis of alcohol intoxication. The accuracy of this device has been established at a Blood Alcohol Concentration of 0.008 and 0.032 % BAC. Accuracy at a Blood Alcohol Concentration greater than 0.032% BAC has not been established.

## Preparation for Correct Measurement


- Before testing, wait at least 15 minutes after your last drinking.
- Only use the device in a temperature 10°C -40°C or 50°F -104°F.

## Description of parts

## Power Supply

- Slide back the battery cover on the back of the device.
- Insert 2 pieces AAA battery, making sure that the “+/-“symbols are correctly.


- If the power is too low, the unit may display a battery symbol  and the device may not use correctly. You need to change a new battery.


## Operation

### (Test mode)

- Insert one of the mouthpieces into the 'blow in' hole.

- Press and hold the TEST button until power on. Release it after power on.

- The WAIT displays on the screen and countdown 10, 9, and 8... to zero.

- The BLOW displays on the screen, and displays "C" with flashing. Take a deep breath, blowing steady into the mouthpiece until beeping stops (about 6 seconds).

- The device displays "WAIT" and flashing the "0" on the screen.

- Your test result will be displayed on the screen after few seconds.

- Your test result will be changed the color if above alarm point.

- It will show error when incorrect blowing. Refer to the chapter of the “Errors and Troubles” in detail. The device will be turned off after indicated error.

**Accessories**

- Mouthpieces
- Manual

## Errors or troubles

Error message or troubles	Description	Remarks
No display	-Press and hold the TEST button until power on. -Check the battery	
Auto Off	-Low battery -Over the period of test stage	
Err	-Blowing unsteady or interrupt -Blowing less than 6 seconds -Blowing flow too weak -No blow during BLOW mode. -Change temperature suddenly	
	-Low battery indication -Need to change new battery -Use alkaline battery	
Alarm sound	-During the blowing mode -Error message indication -Above alcohol set point	
HI	-The blood alcohol content above the display range. -Display range: 0.20%BAC(2.00‰, 1.000mg/L)	

## Important notes

- Before testing, wait at least 15 minutes after your last drinking.
- The blowing time must be enough until beeping stops.
- It may show ERR (error) if none, interrupt or shorter blowing.
- If the device changes temperature suddenly, wait some minutes to use.
- Only use the device in a temperature 10°C -40°C or 50°F -104°F.
- Recalibration interval: 12 months or 200 measurements.
- Never store the device in a corrosive gas room.
- Never store the device under sun exposure area.
- Avoid dropping or subjecting the device to a strong shock.

## Specifications

- Electrochemical sensor (Fuel Cell Sensor)
- Blowing time about 6 seconds
- Blowing flow about 10L/min
- Blowing flow interrupt detection
- Test counter to display the operated number
- 4 digits display with status and EL backlight
- Display range: 0.000-0.200%BAC(2.00‰ or 1.00mg/L)
- Show "HI" if above the display range
- Accuracy: +/-0.005% at 0.05%
- 5 mouthpieces included
- Warm up time: 10-15 seconds
- Testing Time: 10 seconds
- Repeat test: no less than 3 minutes.
- Audible sound and change LCD color if above alarm point
- Power supply: 2 pcs AAA 1.5V alkaline battery
- Battery for standard testing is about 300 times
- Low battery indication and auto power off function
- Operation: 10°C ~ 40°C or 50°F ~ 104°F
- Store: 0°C ~ 50°C or 32°F~122°F
- Size: 115x60x23mm

## Calibration and record

- You should have your alcohol tester calibrated every 12 months or 200 measurements. Failure to do so may cause your readings to be inaccurate.

Calibration Record			
Date First Used		Serial Number	
#	Calibration Date	Remarks	Recalibration by

## Precautions

- Your blood alcohol content (BAC) may change some minutes later.
- You should have your alcohol tester calibrated every 12 months or 200 measurements. Failure to do so may cause your readings to be inaccurate.
- Do not use this device to determine if you can drive.
- You cannot use the test results of device in court of law.
- You cannot use the test result of device against police's unit.

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## English

### **Correct Disposal of This Product (Waste Electrical & Electronic Equipment)**

**(Applicable in the European Union and other European countries with separate collection systems)**

This marking shown on the product or its literature, indicates that it should not be disposed with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

Made in China for LECHPOL Zbigniew Leszek, Miętne, 1 Garwolińska Street 08-400 Garwolin.

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